

Over Toilet Frame



For further information please contact:
06 878 8109 ext 5761



HAWKE'S BAY
District Health Board
Whakawateatia

Reviewed June 2014

Instructions for using an
Over Toilet Frame at
home



HAWKE'S BAY
District Health Board
Whakawateatia

Instructions for using an Over Toilet Frame at home

To adjust height:



- Push in the metal buttons on the legs
- Set legs at the height to suit the person using the frame
- Click the buttons into the correct holes, make sure the buttons are completely popped out
- Make sure all four legs are the same height
- One leg has extra holes—this is used to adjust the frame for uneven surfaces.

To Use:



- Make sure feet are flat on floor and hips are above the knees
- Always use two hands to push off the frames arms as it can tip
- The toilet seat may need to be removed to fit the frame.

DO NOT USE FRAME IN THE SHOWER OR BATH

Maximum weight capacity is 115kg

Additional information from your Occupational Therapist:

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

If you have any questions or require further advice please ring the

Hawke's Bay Hospital Equipment Loan Store
Phone: 06 878 8109 ext 5761